Date 21-09-21		
Expt. No Page No		
DETERMINATION OF Na2CO2 AND NAOH IN A MIXTURE BY TITRATION		
* AIM: To determine the amount of Na2 CO3 and NaOH in a mixture using hydrochloric acid.		
# PRINCIPLE: When a known volume of the mixture is titrated with the in presence of phenolphthalein, the acid reacts with all the sodium hydroxide and with only half of the carbonate. (2) When the titration is continued with methyl wang indicator, the remaining half of CO2 ions will be reutralised with HCl at the end point. A = all hydroxide wins + half of carbonate wins. B = half of carbonate wins after phenolphthalein end point 2B = all carbonate wins. A - B = all hydroxide wins.		
* PROCEDURE: (1) TITRATION 1: Standardization of HCl. -) 20 ml of 0.1 N Na, Co2 solution is pipetted into a clean conical flask and two on three drops of methyl arange indicator is added to the solution. -) The selution is the titrated against hydrochloric acid taken in the burette. -> The end point is change of colour from yellow to crange. The titrations are repeated to get concordant values		

	No.	
Expt. NoPage No		
10)	TITRATION D. O.A. L.	
7	THE given unknown addition of Na, CO3 and NaOH.	
	The given unknown solution is made upto 100 ml	
2	a standard flask using distilled water.	
->	and of this made up solution is bibetted into	
_	in a standard plack using distilled water. 20 ml of this made up solution is pipetted into a clean conical plack. Add 2 to 3 drops of phenolphalein indicator and titrate against standard MCL.	
-		
	titrate against standard HCL.	
->	The ansappearance of bink colour is also	
	foint and note down the titre value. Let it be A ml.	
-	solution add 2 or 3 dide I will	
	changes from yellow to orange and note down the titre value after the phenolphthalein and point. Let it be B ml.	
	value after the phenolphtlalein and boint. Let it be	
	Bril.	
->	The titration is repeated to get concordant values of both A and B.	
	both A and B.	
*	OBSERVATIONS AND CALCULATIONS:	
7	KNOWN VALUES:	
	Volume of HCl (V,) = 24.4 ml	
	Normality of MCl (N1) = ?	
	volume of Na, CO3 (V2) = 20 ml	
	Normality of Na, co3 (N2) = 0.1 N	
	Volume of Kel required for neutralization of Naz CO3 = 18 ml	
	:. Normality of HCl, N, = 20 x 0.1 = 20 x 0.1 > 0.081967 N	
	V ₁ 24.4	
	: Normality of Naylo3 = 18 x 0.08 1967 => 0.07377 N	
	20 Teacher's Signature	