

Expt. - Determination of molecular weight of a polymer by viscosity Average method

Aim

To determine the molecular weight of a polymer in solution by using a viscometer.

Apparatus Required

Ostwald's viscometer, volumetric flasks stop watch, standard flasks.

Reagent required

Polymer, suitable solvents.

Principle

If a polymer is soluble in a suitable solvent, measurement of solution viscosity provides a simple and convenient method. For a particular molecular weight determination. In a capillary viscometer (Ostwald/Ubbelohde) the viscosity of a liquid is proportional to the time taken by a known volume of a liquid to flow through a capillary under a specified hydrostatic pressure at a fixed

temperature. The condition for the flow should ensure that the flow is laminar. Using Poiseuille's equation it is possible to show that t, η and ρ are the flow time, viscosity and density of a solution respectively; and t_0, η_0, ρ_0 are those of the pure solvent, then

$$\eta/\eta_0 = \rho/\rho_0 \cdot t/t_0$$

The value of η/η_0 , is known as the relative viscosity, η_{rel} . In dilute solutions, which are often employed. For molecular weight determination, ρ is not much different from ρ_0 and hence

$$\eta_{rel} = \eta/\eta_0 = t/t_0$$


The specific viscosity is defined as

$$\eta_{sp} = \eta_{rel} - 1$$

A plot of η_{sp}/c vs c is a straight line for dilute solutions, the intercept

$$c \xrightarrow{\text{lim}} 0 \left(\frac{\eta_{sp}}{c} \right) = \eta_{int}$$

Aim → To determine the molecular weight of a polymer in solution by using a viscometer



- Plot of (η_{sp}/c) vs Concentration of Polymer solution to find out intrinsic viscosity.

* Calculate η_{sp} and η_{sp}/c . Plot η_{sp}/c vs c , extrapolate to $c=0$ to obtain η_{int} . From the given values of K and α , calculate the molecular weight.

Sr. No.	Concentration of the polymer solution	Time of flow in sec (average)	Relative viscosity $\eta_{rel} = \frac{\eta}{\eta_0} = t/t_0$	Specific viscosity $\eta_{sp} = \eta_{rel} - 1$	Reduced viscosity η_{sp}/c
1.	Pure Solvent	$t_0 = 55$	—	—	—
2.	0.1%	$t_s = 58$	1.0545	0.0545	54.5
3.	0.2%	$t_s = 62$	1.1273	0.1273	63.6
4.	0.3%	$t_s = 67$	1.2181	0.2181	72.7
5.	0.4%	$t_s = 72$	1.3090	0.3090	77.2
6.	0.5%	$t_s = 79$	1.4363	0.4363	87.2

$$[K = 45.3 \times 10^{-3}, \alpha = 0.64]$$


Charvi Jain

RA2111047010113

Scale

X-axis : 2 cm - 0.1%

Y-axis : 2 cm - 10 units



Concentration of Polymer solution
(C%) \longrightarrow

* Calculation

solvent used = water

$$\text{Now, } \eta_{\text{int}} = k \times M^{\alpha}$$

Taking log on both sides

$$\log \eta_{\text{int}} = \log k + \log M^{\alpha}$$

$$\Rightarrow \log \eta_{\text{int}} = \log k + \alpha \log M$$

$$\Rightarrow \alpha \log M = \log \eta_{\text{int}} - \log k$$

$$\text{From, } \eta_{\text{int}} = 46.5, \text{ given } k = 45.3 \times 10^{-9}, \\ \alpha = 0.64$$

$$\Rightarrow \log M = \frac{\log 46.5 - \log 45.3 \times 10^{-9}}{\alpha}$$

$$= \frac{\log 46.5 - \log 45.3 + 3 \log 10}{0.64}$$

$$= \frac{1.66 - 1.65 + 3}{0.64} = \frac{3.01}{0.64} = 4.703$$

$$\therefore M = \text{antilog}(4.703) = 50,466.13 \text{ g}$$

* Result

The molecular weight of given polymer = 50,466.13 g

which is known as the intrinsic viscosity η_{int} .

The Staudinger - Mark - Houwink equation which relates η_{int} with molecular weight

$$\eta_{int} = k(M)^{\alpha}$$

where 'k' is an empirical parameter characteristic of a particular solute-solvent pair and 'α' is a 'shape' parameter, characteristic which can vary from about 0.5 for well coiled polymers in poor solvents to about 2. For rigidity extended 'rod' like polymers. From known values of k and α, Molecular weight can be determined.

Procedure

Preparation of variation concentration of polymer in water (solvent)

1% solution of polymer in water will be supplied. We need to prepare at least '5' dilution viz. 0.1%, 0.2%, 0.3%, 0.4% and 0.5% polymer in water before carrying out the experiment.

Dilutions can be done by using volumetric solution from a 1% solution, volume is

$$V_1 = \frac{V_2 \times N_2}{N_1} = \frac{100 \text{ ml} \times 0.2\%}{1\%} = 20 \text{ ml}$$

Similarly, any other dilution can be prepared by the above method.

Setup the Ostwald (or Ubbelohde) viscometer and measure the flow time (t_0) of a fixed volume of the pure solvent. Take an average of three readings. Rinse the viscometer thoroughly with the most dilute solution, measure the flow time (t_1) keeping the flow-volume the same. Repeat the procedure for other solutions.