


(1) Pilot titration



(2) Fair titration

## ESTIMATION OF IRON BY POTENTIOMETRIC TITRATION

### \* AIM:

To estimate the amount of  $\text{Fe}^{2+}$  ion present in the given solution.

### \* APPARATUS REQUIRED:

Potentiometer assembly, 25 ml burette, 10 ml pipette, 250 ml beakers, standard flask, calomel and platinum electrodes.


### \* REAGENTS REQUIRED:

Ferrous ammonium sulphate, dil  $\text{H}_2\text{SO}_4$ , std.  $\text{K}_2\text{Cr}_2\text{O}_7$ .

### \* PRINCIPLE:

→ In a titration, there is a change in ionic concentration which can be followed by measuring the potential of a suitable electrode.

→ The potential of an electrode depends upon the concentration of the ions in accordance to Nernst equation,




$$\therefore E = E^\circ - (RT/nF) \ln [\text{R}/\text{O}^{n+}]$$

where  $E^\circ$  is standard reduction potential.

→ The cell setup:  $\text{Hg}|\text{HgCl}_2(\text{s}), \text{KCl}(1\text{N})||\text{Fe}^{2+}|\text{Fe}^{3+}, \text{Pt}$

→ The chemical reaction:



→

P.T.O

Teacher's Signature \_\_\_\_\_

★ OBSERVATIONS:

(1) TABLE 1:

S.No	Vol. of $K_2Cr_2O_7$ (ml)	EMF (V)	$\Delta E$ (V)	$\Delta V$ (ml)	$\Delta E / \Delta V$ ( $V_{ml}^{-1}$ )
1	7.0	433	—	—	—
2	7.2	441	8	0.2	40
3	7.4	452	11	0.2	55
4	7.6	464	12	0.2	60
5	7.8	516	52	0.2	260
6	8.0	599	83	0.2	415
7	8.2	622	23	0.2	115
8	8.4	632	10	0.2	50
9	8.6	640	8	0.2	40
10	8.8	645	5	0.2	25
11	9.0	649	4	0.2	20
12	9.2	652	3	0.2	15
13	9.4	655	3	0.2	15
14	9.6	658	3	0.2	15
15	9.8	660	2	0.2	10
16	10.0	667	2	0.2	10



### \* PROCEDURES:

- The given ferrous ion solution is made up to a known volume in a SMF following the standard procedure with usual precautions.
- Exactly 10 ml of the made up  $\text{Fe}^{2+}$  solution is pipette out into a clean 100 ml beaker. About 10 ml of dilute  $\text{H}_2\text{SO}_4$  and 100 ml of distilled water are added to it.
- A platinum electrode is dipped into this solution and coupled with the standard calomel electrode. The resultant cell is then incorporated into the potentiometric circuit.
- Standard  $\text{K}_2\text{Cr}_2\text{O}_7$  solution, which is taken in a burette, is added in installments of 1 ml into the beaker and the cell e.m.f is measured after each addition by proper mixing.
- The process is continued till and also well beyond the neutralization point as indicated by an abrupt change in the emf.
- Note the volume  $\text{K}_2\text{Cr}_2\text{O}_7$  solution required for complete oxidation of  $\text{Fe}^{2+}$  solution from the plot of e.m.f versus the volume of the standard  $\text{K}_2\text{Cr}_2\text{O}_7$  solution added. The range at which the endpoint lies maybe evaluated.
- Calculate the normality of given  $\text{Fe}^{2+}$  solution using the formula  $N_1V_1 = N_2V_2$
- One more similar titration is performed by adding 0.1 ml portions of standard  $\text{K}_2\text{Cr}_2\text{O}_7$  solution close to the end point and tabulating the measured e.m.f corresponding to each addition.

Teacher's Signature \_\_\_\_\_

(2) TABLE 2:

S.No.	Vol. of $K_2Cr_2O_7$ (ml)	EMF (V)	$\Delta E$ (V)	$\Delta V$ (ml)	$\Delta E/\Delta V$ (Vml <sup>-1</sup> )
1	7.0	433	—	—	—
2	7.2	441	8	0.2	40
3	7.4	452	11	0.2	55
4	7.6	464	12	0.2	60
5	7.8	576	52	0.2	260
6	8.0	599	83	0.2	415
7	8.2	622	23	0.2	115
8	8.4	632	10	0.2	50
9	8.6	640	8	0.2	40
10	8.8	645	5	0.2	25
11	9.0	649	4	0.2	20
12	9.2	652	3	0.2	15
13	9.4	655	3	0.2	15
14	9.6	658	3	0.2	15
15	9.8	660	2	0.2	10
16	10.0	667	2	0.2	10



- Graph is plotted (1) EMF vs Volume of  $K_2Cr_2O_7$  addition.  
(2)  $\Delta E/\Delta V$  vs Volume of  $K_2Cr_2O_7$  addition.
- The exact end point can be determined from the plot of  $\Delta E/\Delta V$  vs Volume of standard  $K_2Cr_2O_7$  solution.

### ★ CALCULATIONS:

→ Volume of pipette solution (FAS) =  $V_1$  ml = 10 ml

Volume of  $K_2Cr_2O_7$  =  $V_2$  ml (from graph) = 8 ml

Normality of  $K_2Cr_2O_7$  =  $N_2$  = 0.1 N

Normality of FAS ( $N_1$ ) = ?

$$\text{Strength of FAS} = N_1 = \frac{N_2 V_2}{V_1} = \frac{8 \times 0.1}{10} = \frac{0.8}{10} = 0.08 = 0.008 N //$$

$$\text{Amount of } Fe^{2+} (g\ l^{-1}) = 55.85 \times 0.008 \\ = 0.4468$$

$$\therefore \text{Amount of } Fe^{2+} \text{ in } 100\text{ ml} = \frac{0.4468}{10} \Rightarrow 0.04468\text{ g} //$$




P.T.O

Teacher's Signature \_\_\_\_\_



# GRAPHS:






# GRAPHS:

Expt. No.

★ RE

→ TU





★ RESULT:

→ The weight of Fe present in the given solution = 0.4468g

— X —

Teacher's Signature \_\_\_\_\_