

Experiment-2 : Conductometric titration

Conductance measurements are frequently used to find the end points of acid alkali and other titrations.

Principle: The electrical conductance depends upon the number and mobility of ions

Dt. 26/4/2021


Conductometric titration

Case-1: A strong acid against a strong base

Ex: HCl vs. NaOH



- ❖ On addition of NaOH, the H^+ ions are replaced by slow moving Na^+ ions.
- ❖ After neutralization, addition of alkali will introduce fast moving hydroxyl ions.




Case-2: A weak acid against a strong base

Ex: CH_3COOH vs. NaOH




On adding the base, highly ionized Sodium acetate is formed and hence the conductance begins to increase

After neutralization, addition of the base will introduce fast moving hydroxyl ions.



Case-3: A mixture of strong and weak acid against a strong base

Ex: Mixture of HCl and CH_3COOH vs. NaOH




Case-4: A strong acid against a weak base

Ex: **HCl** vs. **NH₄OH**



- ❖ The conductance will fall first on addition of NaOH, because the fast moving H^+ ions are replaced by slow moving Na^+ ions.
- ❖ After neutralization of the acid, further addition of weakly ionized NH_4OH will not cause any appreciable change in the conductance.



Advantages

- ❖ Coloured solutions which can't be titrated by ordinary volumetric methods with the help of Indicators, can be successfully titrated conductometrically
- ❖ This method can be employed in case of very dilute solutions and also for weak acids and bases