Experiment No. 1

Aim of the experiment:

To determine the amount of Na₂CO₃ and NaOH in a mixture using hydrochloric acid.

Principle:

When a known volume of the mixture of Na₂CO₃ and NaOH is titrated with HCl using phenolphthalein indicator, at the end point, all the hydroxide ions (OH⁻) and only half of the carbonate ions (CO₃²⁻) are reacted with hydrochloric acid.

When the titration is continued with methyl orange indicator, the remaining half of CO_3^{2-} ions will be neutralized with HCl at the end point

$$NaOH + HC1$$
 — $NaCl + H2O$

$$Na_2CO_3 + 2 HCl$$
 $2NaCl + H_2O + CO_2$

1st Titration

Titration of HCl against Na₂CO₃

Aim: To know the strength of HCl solution

Performing chemical reactions quantitatively to determine the exact amount of a reagent is called a **titration**

A titration is the quantitative reaction of an acid and a base.

Indicators are used to show that all the analyte has reacted with the titrant.

Known things:

Concentration of Na_2CO_3 solution: $N_1 = 0.1$

Volume of Na₂CO₃ solution: V₁ (20 ml)

Known volume of HCl consumed (from titration): V₂

Concentration of HCl solution: N_2 ??

Strength of HCl solution,

$$N_2 = N_1 V_1 / V_2$$

= 0.1x20/24.4 = 0.081967 N

2nd Titration

Titration of HCl against mixture of NaOH and Na₂CO₃

A = volume of HCl required for neutralization of hydroxyl as well as half of carbonate ions

B= volume of HCl required for neutralization of whole carbonate ions

B-A = C = neutralization of half carbonate ions Volume required for neutralization of whole carbonate ions= 2C ml

Volume of HCl required for neutralization of only hydroxyl ions = A-C m

Calculations

Known things:

Volume of HCl required for neutralization of Na₂CO₃=2 C = 2 (B-A) ml = 2(35.3-26.3) =18 ml

Strength of HCl = 0.081967 N

Volume of $Na_2CO_3 = 20 \text{ ml}$

Strength of $Na_2CO_3 = 18x0.081967/20$ =0.07377 N

Amount of Na_2CO_3 present in the mixture = 0.07377x53/10 = 0.39098 g

Volume of HCl required for C = B - A = 9 ml

Volume of HCl required for neutralization of only hydroxyl ions,

$$A-C = 26.3-9 = 17.3 \text{ ml}$$

Strength of hydroxyl ions = 0.081967x 17.3/20

= 0.0709 N

Amount of NaOH present in the mixture = $0.0709 \times 40/10$

= 0.2836 g

Conclusion

The amount of Na₂CO₃ and NaOH present in the mixture are found to be 0.39098 g and 0.2836 g, respectively.

Normality of NaOH = x

X gm equivalent of NaOH is present in 1 litre or 1000 ml of the solution

1 gm equivalent of NaOH = 40

X gm equivalent = 40 . X gm of NaOH is present in 1000 ml of the solution

40.x/10 gm of NaOH is present in 100 ml of the solution