

1. A mixture of 0.1015 g of NaCl and 0.1324 g of KCl is dissolved in water and titrated with 0.15 M of AgNO<sub>3</sub> (aq.).

What volume of the AgNO<sub>3</sub> (aq.) will be needed to reach the endpoint?

Ans: Milli gm equivalent of NaCl =  $(0.1015/58.46) \times 1000 = 1.736$

Milli gm equivalent of KCl =  $(0.1324/74.5) \times 1000 = 1.777$

Total milli gm equivalent of the mixture =  $1.736 + 1.777 = 3.513$

Concentration of AgNO<sub>3</sub> = 0.15 M

At equivalent point, total milli gm equivalent of chloride = total milli gm equivalent of AgNO<sub>3</sub>

Milli gm equivalent of AgNO<sub>3</sub> = conc. of AgNO<sub>3</sub> x volume of AgNO<sub>3</sub>

$$\text{Hence } 3.513 = 0.15 \times (x \text{ ml})$$

$$\text{or, } x = 3.513/0.15 = 23.42 \text{ ml}$$

Normality = gm equivalent of that substance/ volume

2. How many ml of 0.45 M HCl must be added to 25 ml of 1.0 M KOH to make a neutral solution?

3. What is the cell potential for a Zn/Cu cell when  $[\text{Zn}^{2+}] = 10 \text{ M}$  and  $[\text{Cu}^{2+}] = 1 \text{ M}$  at  $25^\circ\text{C}$ , where for  $\text{Cu}^{2+}(\text{aq}) + 2\text{e}^- \rightarrow \text{Cu}(\text{s})$ ,  $E^\circ = + 0.34 \text{ V}$  and  $\text{Zn}(\text{s}) \rightarrow \text{Zn}^{2+}(\text{aq}) + 2\text{e}^-$   $E^\circ = + 0.76 \text{ V}$ .



$$E(\text{Zn}^{2+}/\text{Zn}) = E^\circ - (0.059/2) \log 1/10$$

$$= -0.76 - 0.0295 (-1)$$

$$= -0.76 + 0.0295$$

$$= -0.7305 \text{ V}$$

$$E_{\text{cell}} = E(\text{Cu}^{2+}/\text{Cu}) - E(\text{Zn}^{2+}/\text{Zn})$$

$$= 0.34 - (-0.7305)$$

$$= 0.34 + 0.7305$$

$$= 1.0705 \text{ V}$$

4. 25 ml of 0.01 M AgNO<sub>3</sub> solution is mixed with 25 ml of 0.0005 M aqueous NaCl solution. Whether there will be any precipitate? Given K<sub>sp</sub> of AgCl is  $1.7 \times 10^{-10} \text{ M}^2$

Ans: [Ag<sup>+</sup>] = 0.005 M

[Cl<sup>-</sup>] = 0.00025 M

Ionic product of AgCl = [Ag] [Cl<sup>-</sup>]

$$= 5 \times 10^{-3} \times 2.5 \times 10^{-4}$$

$$= 12.5 \times 10^{-7}$$

$$= 1.25 \times 10^{-6} \text{ M}^2$$

As the ionic product exceeds the solubility product, there will be precipitation of AgCl.

**5. 10 ml of 0.1 N NaOH is added to 20 ml of 0.1N H<sub>2</sub>SO<sub>4</sub> and the resultant solution is titrated against 0.1 N NaOH solution. What will be the titre value at the end point?**

**Ans: milli gm Equivalent of NaOH =  $0.1 \times 10 = 1$**

**milli gm equivalent of H<sub>2</sub>SO<sub>4</sub> =  $20 \times 0.1 = 2$**

**Net milli gm equivalent of H<sub>2</sub>SO<sub>4</sub> =  $2 - 1 = 1$**

**The net milli gm equivalent of NaOH should be 1**

**Normality of NaOH x volume of NaOH (x) = 1**

**Thus,  $0.1 \times x = 1$**

**Or,  $x = 1/0.1 = 10$  ml**

6. The molar conductance of  $\text{CH}_3\text{COONa}$ ,  $\text{HCl}$  and sodium chloride at infinite dilution are  $91 \times 10^{-4}$ ,  $426.16 \times 10^{-4}$  and  $126.45 \times 10^{-4} \text{ S m}^2 \text{ mol}^{-1}$ , respectively at  $25^\circ\text{C}$ . What will be the molar conductance value at infinite dilution for acetic acid?

Ans:  $\text{CH}_3\text{COONa} + \text{HCl} \longrightarrow \text{CH}_3\text{COOH} + \text{NaCl}$

$$\text{CH}_3\text{COOH} = \text{CH}_3\text{COONa} + \text{HCl} - \text{NaCl}$$

$$= (91 + 426.16 - 126.45) \times 10^{-4}$$

$$= 390.71 \times 10^{-4} \text{ S m}^2 \text{ mol}^{-1}.$$

7. In determination of mixture of bases by titration method, the amount of Sodium Hydroxide is calculated as---.

- a.  $N \times \text{Equivalent mass of Sodium Carbonate} / 10$
- b.  $N [\text{OH and CO}_3^{2-} \text{ portion}] \times \text{Equivalent mass of Sodium Hydroxide and Sodium carbonate} / 10$
- c.  $N [\text{OH portion}] \times \text{Equivalent mass of Sodium Hydroxide} / 10$
- d.  $N [\text{CO}_3^{2-} \text{ portion}] \times \text{Equivalent mass of Sodium carbonate} / 10$

8. What is the working principle of conductometry?

- a. measurement of potential.
- b. measurement of conductivity of solution.
- c. measurement of emf.
- d. measurements of pH

8. Conductivity cell is made up of...

- a. Two silver rods
- b. Two parallel sheets of platinum
- c. Glass membrane of Ag/AgCl
- d. Sb-Sb<sub>2</sub>O<sub>3</sub>

9. At the same concentration and temperature, dilute aqueous solution of strong acid will conduct electricity....

- a. better than dilute aqueous solution of weak acid
- b. as much as dilute aqueous solution of weak acid
- c. lower than the dilute aqueous solution of weak acid
- d. two-fold higher than the weak acid



10. If 2 g of NaOH is dissolved in 100 ml of distilled water, then what is the concentration of the solution?

a. 1 N

b. 2 N

c. 0.5 N

d. 0.05 N

11. Which of the following is the formula for pH calculation?

a)  $\log_{10}[\text{H}^+]$

b)  $-\log_{10}[\text{H}^+]$

c)  $\log_2[\text{H}^+]$

d)  $-\log_2[\text{H}^+]$

12. The pH meter is a

- a) Ammeter
- b) Voltmeter
- c) Potentiostat
- d) Spectrophotometer

**13. How we will come to know that a given solution is acidic?**

- a) If its pH value is less than 7
- b) If its pH value is greater than 7
- c) If its pH value is less than 5
- d) If its pH value is 5

**14. A buffer solution is used with pH measuring instruments to**

- a) protect the equipment
- b) standardize the equipment
- c) clean the electrodes
- d) platinize the reference electrode

**15. The pH of a liquid solution is a measure of**

- a) dissolved salt content
- b) hydrogen ion activity
- c) hydroxyl ion molarity
- d) electrical conductivity

16. Measurement of solution viscosity offers a simple and convenient method for molecular weight determination if

- a) Polymer is insoluble in solvent
- b) Polymer is soluble in solvent
- c) Polymer is sparingly soluble in solvent
- d) Polymer is used as neat

17. Which one of the following equations is used to calculate the relative viscosity?

- a)  $\eta / \eta_0 = t/t_0$
- b)  $\eta_{sp} = \eta / \eta_0 - 1$
- c)  $\eta_{red} = \eta_{sp} / C \times 100$
- d)  $\eta_i = K (M)^a$

18. Viscosity is due to one of the following

- a) Potential energy stored in fluid
- b) Resistance to fluid motion
- c) Roughness of the surface
- d) The pressure difference between the two fluids

19. What is the role of chromate ions in chloride estimation?

- a. It acts as a reducing agent
- b. It acts as a buffer
- c. It acts as an indicator
- d. It acts as an oxidizing agent

20. What is the pH range in which chloride determination using Mohr's method is conducted?

- a.  $< 3$
- b. 5
- c.  $> 12$
- (d) 6 -9

21. Why do we have to standardize  $\text{AgNO}_3$  solution?

- a. To find the normality of  $\text{NaCl}$
- b. To calculate the normality of  $\text{AgCl}$
- c. To find the normality of  $\text{AgNO}_3$
- d. To calculate the volume of  $\text{NaCl}$

22. All of the following statements are correct regarding potentiometric titration except

- a. They are suitable for colored or turbid solutions
- b. The EMF of the cell is zero at the equivalence point
- c. The results obtained are accurate
- d. Acid base titration can also be carried out by potentiometry

23. Hard water + Buffer + EBT -----□

- a. Appearance of wine-red colour
- b. Appearance of steel blue colour
- c. Formation of weak complex
- d. Formation of brown precipitate

24. Hard water + Buffer + EBT + EDTA

- a. Appearance of wine-red colour
- b. Appearance of steel blue colour
- c. Formation of weak complex
- d. Formation of brown precipitate